

Anti-MTHFD2L antibody (STJA0004787)

STJA0004787

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Bifunctional Methylenetetrahydrofolate Dehydrogenase/Cyclohydrolase 2-Mitochondrial is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.84 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:30-1:150
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	441024
Gene Symbol	MTHFD2L
Uniprot ID	MTD2L_HUMAN
Immunogen	Fusion protein of human MTHFD2L
Immunogen Region	
Specificity	
Immunogen Sequence	