

Anti-KRCC1 antibody (STJA0004808)

STJA0004808

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Lysine-Rich Coiled-Coil Protein 1 is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
Applications	WB/IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.68 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	WB 1:500-1:2000, IHC 1:100-1:300
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51315
Gene Symbol	KRCC1
Uniprot ID	KRCC1_HUMAN
Immunogen	Fusion protein of human KRCC1
Immunogen Region	
Specificity	
Immunogen Sequence	