

Anti-C1orf112 antibody (STJA0004823)

STJA0004823

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Uncharacterized Protein C1orf112 is suitable for use in Immunohistochemistry and ELISA research applications. |
| Applications | IHC/ELISA |
| Host/Source | Rabbit |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.96 mg/mL |
| Conjugation | Unconjugated |
| Purification | Antigen affinity purification |
| Dilution Range | IHC 1:50-1:100 |
| Formulation | PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4 |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term, and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|----------------------------------|
| Gene ID | 55732 |
| Gene Symbol | FIRRM |
| Uniprot ID | FIRRM_HUMAN |
| Immunogen | Fusion protein of human C1orf112 |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |