

Anti-EXD2 antibody (STJA0004825)

STJA0004825

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Rabbit polyclonal antibody anti-Exonuclease 3 5 Domain-Containing Protein 2 is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications. |
| Applications | WB/IHC/ELISA |
| Host/Source | Rabbit |
| Reactivity | Human/Mouse |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 0.96 mg/mL |
| Conjugation | Unconjugated |
| Purification | Antigen affinity purification |
| Dilution Range | WB 1:500-1:2000, IHC 1:40-1:200 |
| Formulation | PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4 |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for long term, and avoid freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|------------------------------|
| Gene ID | 55218 |
| Gene Symbol | EXD2 |
| Uniprot ID | EXD2_HUMAN |
| Immunogen | Fusion protein of human EXD2 |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |