

Anti-GOLGA6A antibody (STJA0004835)

STJA0004835

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Golgin Subfamily A Member 6a is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.38 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:50-1:300
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	342096
Gene Symbol	GOLGA6A
Uniprot ID	GOG6A_HUMAN
Immunogen	Synthetic peptide of human GOLGA6A
Immunogen Region	
Specificity	
Immunogen Sequence	