

Anti-CAMSAP3 antibody (STJA0004847)

STJA0004847

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Calmodulin-Regulated Spectrin-Associated Protein 3 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1.32 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	IHC 1:50-1:200
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	57662
Gene Symbol	CAMSAP3
Uniprot ID	CAMP3_HUMAN
Immunogen	Synthetic peptide of human CAMSAP3
Immunogen Region	
Specificity	
Immunogen Sequence	