

Anti-PEMT antibody (STJA0004854)

STJA0004854

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phosphatidylethanolamine N-Methyltransferase is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	2.4 mg/mL
Conjugation	Unconjugated
Purification	Antigen affinity purification
Dilution Range	WB 1:1000-1:5000
Formulation	PBS with 0.05% NaN ₃ and 40% Glycerol, pH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10400
Gene Symbol	PEMT
Uniprot ID	PEMT_HUMAN
Immunogen	Synthetic peptide of human PEMT
Immunogen Region	
Specificity	
Immunogen Sequence	