

Anti-TRIM10 antibody (STJA0004855)

STJA0004855

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Tripartite Motif-Containing Protein 10 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:25-100
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	10107
Gene Symbol	TRIM10
Uniprot ID	TRI10_HUMAN
Immunogen	Fusion protein of human TRIM10
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	