

**Anti-DUSP13 antibody (STJA0004859)**

STJA0004859

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-Dual Specificity Protein Phosphatase 13 Isoform B is suitable for use in Immunohistochemistry and ELISA research applications.
<b>Applications</b>	IHC/ELISA
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	0.3 mg/mL
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	IHC 1:50-1:200
<b>Formulation</b>	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for long term, and avoid freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	51207
<b>Gene Symbol</b>	DUSP13B
<b>Uniprot ID</b>	DS13B_HUMAN
<b>Immunogen</b>	Recombinant protein of human DUSP13
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	