

Anti-ENPP5 antibody (STJA0004864)

STJA0004864

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ectonucleotide Pyrophosphatase/Phosphodiesterase Family Member 5 is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.5 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:50-1:200
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	59084
Gene Symbol	ENPP5
Uniprot ID	ENPP5_HUMAN
Immunogen	Recombinant protein of human ENPP5
Immunogen Region	
Specificity	
Immunogen Sequence	