

Anti-SSX2IP antibody (STJA0004881)

STJA0004881

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Afadin-And Alpha-Actinin-Binding Protein is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
Applications	WB/IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.7 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000, IHC 1:40-1:150
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	117178
Gene Symbol	SSX2IP
Uniprot ID	ADIP_HUMAN
Immunogen	Recombinant protein of human SSX2IP
Immunogen Region	
Specificity	
Immunogen Sequence	