

Anti-SSX2IP antibody (STJA0004882)

STJA0004882

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Afadin-And Alpha-Actinin-Binding Protein is suitable for use in Immunohistochemistry and ELISA research applications.
Applications	IHC/ELISA
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.6 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC 1:30-1:150
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage Instruction	Store at -20°C for long term, and avoid freeze-thaw cycles.

TARGET INFORMATION

Gene ID	117178
Gene Symbol	SSX2IP
Uniprot ID	ADIP_HUMAN
Immunogen	Recombinant protein of human SSX2IP
Immunogen Region	
Specificity	
Immunogen Sequence	