

Anti-AUP1 antibody (STJA0004896)

STJA0004896

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Lipid Droplet-Regulating Vldl Assembly Factor Aup1 is suitable for use in Western Blot, Immunohistochemistry and ELISA research applications.
Applications	WB/IHC/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	0.4 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000, IHC 1:50-1:200
Formulation	PBS with 0.05% sodium azide and 50% glycerol, PH7.4
Isotype	IgG
Storage	Store at -20°C for long term, and avoid freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	550
Gene Symbol	AUP1
Uniprot ID	AUP1_HUMAN
Immunogen	Recombinant protein of human AUP1
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	