

Anti-CNP antibody (Full Length) (STJA0004898)

STJA0004898

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Goat polyclonal antibody anti-2, 3-cyclic nucleotide-3-phosphodiesterase (Full Length) is suitable for use in Western Blot and Immunohistochemistry research applications.
Applications	WB/IHC
Host/Source	Goat
Reactivity	Bovine/Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Affinity Purified
Dilution Range	WB 1:5000
	IHC 1:3000
Formulation	PBS + 50% glycerol and 5 mM Na ₃
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1267
Gene Symbol	CNP
Uniprot ID	CN37_HUMAN
Immunogen	Full length human recombinant protein
Immunogen Region	Full Length
Specificity	Specific for endogenous levels of the ~46 kDa CNP protein.
Immunogen Sequence	Full-Length