

Anti-MeCP2 antibody (Full Length) (STJA0004901)

STJA0004901

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Chicken polyclonal antibody anti-MeCP2 (Full Length) is suitable for use in Western Blot and Immunohistochemistry research applications.
Applications	WB/IHC
Host/Source	Chicken
Reactivity	Human/Mouse/Non-Human Primates/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Total IgY fraction
Dilution	WB 1:20000
Range	IHC 1:1000-1:2000
Formulation	PBS + 5 mM NaN3.
Isotype	IgY
Storage	Recommended that the undiluted antibody be aliquoted into smaller working volumes (10-30 uL/vial depending on usage) upon arrival
Instruction	and stored long term at -20°C C or -80°C C, while keeping a working aliquot stored at 4°C C for short term. Avoid freeze/t

TARGET INFORMATION

Gene ID	4204
Gene Symbol	MECP2
Uniprot ID	MECP2_HUMAN
Immunogen	Full length human recombinant MeCP2
Immunogen Region	Full Length
Specificity	Specific for endogenous levels of the MeCP2 protein. Due to usually high charge, the protein runs at ~75 kDa in SDS-PAGE of mouse brain lysates.
Immunogen Sequence	Full-Length