

Anti-Prosurfactant Protein C antibody (N-Term) (STJA0004903)

STJA0004903

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Prosurfactant Protein C (N-Term) is suitable for use in Immunohistochemistry research applications.
Applications	IHC
Host/Source	Rabbit
Reactivity	Mouse/Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	IHC 1:100
Formulation	PBS
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6440
Gene Symbol	SFTPC
Uniprot ID	PSPC_HUMAN
Immunogen	Synthetic peptide from the N-terminal region of mouse Prosurfactant Protein C.
Immunogen Region	N-Term
Specificity	Specific for endogenous levels of the ~21 kDa Prosurfactant Protein C.
Immunogen Sequence	N-Term