

**Anti-SARM1 antibody (Full Length) (STJA0004904)**

STJA0004904

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-SARM1 (Full Length) is suitable for use in Western Blot and Immunocytochemistry research applications.
<b>Applications</b>	WB/ICC
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human/Mouse/Rat

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Antigen Affinity Purified from Pooled Serum
<b>Dilution Range</b>	WB 1:1000 ICC 1:50
<b>Formulation</b>	10 mM HEPES (pH 7.5) , 150 mM NaCl, 100 µg per ml BSA and 50% glycerol.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	23098
<b>Gene Symbol</b>	SARM1
<b>Uniprot ID</b>	SARM1_HUMAN
<b>Immunogen</b>	Full length recombinant human SARM1 protein
<b>Immunogen Region</b>	Full Length
<b>Specificity</b>	Specific for endogenous levels of the ~79 kDa SARM1 protein.
<b>Immunogen Sequence</b>	Full-Length