

Anti-CD105 antibody [SM4908] (STJA0004908)

STJA0004908

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD105 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM4908
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (M) 1:500-1:1000
Formulation	PBS with 0.02% sodium azide, 100 µg/ml BSA and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	13805
Gene Symbol	Eng
Uniprot ID	EGLN_MOUSE
Immunogen	Recombinant protein corresponding to Mouse CD105
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081