

## Anti-eNOS antibody [SM4915] (STJA0004915)

STJA0004915

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-eNOS is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
<b>Applications</b>	IHC/IF
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Rat

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	SM4915
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	IHC/IF (H, R) 1:1000-1:2000/1:500-1:1000
<b>Formulation</b>	PBS with 0.02% sodium azide, 100 µg/ml BSA and 50% glycerol.
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to one year, and avoid repeated freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	18127
<b>Gene Symbol</b>	Nos3
<b>Uniprot ID</b>	NOS3_MOUSE
<b>Immunogen</b>	KLH conjugated Synthetic peptide corresponding to Mouse eNOS
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	



Immunohistochemistry analysis of paraffin-embedded human liver cancer using eNOS at dilution of 1:1000

Immunohistochemistry analysis of paraffin-embedded human colon cancer using eNOS at dilution of 1:1000

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081