

Anti-Caspase-9 antibody [SM4940] (STJA0004940)

STJA0004940

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Caspase-9 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Mouse


PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM4940
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB (M) 1:100-1:200
Formulation	PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol.
Isotype	IgG2a
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12371
Gene Symbol	Casp9
Uniprot ID	CASP9_MOUSE
Immunogen	KLH conjugated Synthetic peptide corresponding to Mouse Caspase-9
Immunogen Region	
Specificity	
Immunogen Sequence	

Mouse skeletal muscle



Western blot analysis of Caspase-9 at dilution of 1:200

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081