

Anti-smooth muscle Myosin heavy chain 11 antibody [SM5002] (STJA0005002)

STJA0005002

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-smooth muscle Myosin heavy chain 11 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM5002
Concentration	
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution Range	IHC/IF (H) 1:200-1:500
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	
Storage	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	4629
Gene Symbol	MYH11
Uniprot ID	MYH11_HUMAN
Immunogen	Synthesized peptide derived from human MyHC
Immunogen Region	
Specificity	
Immunogen Sequence	

