

Anti-delta 1 Catenin/CAS antibody [SM5005] (STJA0005005)

STJA0005005

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-delta 1 Catenin/CAS is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM5005
Concentration	
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution Range	IHC/IF (H) 1:200-1:500
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	
Storage	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1500
Gene Symbol	CTNND1
Uniprot ID	CTND1_HUMAN
Immunogen	Synthesized peptide derived from human p120 catenin
Immunogen	
Region	
Specificity	
Immunogen	
Sequence	