

Anti-S100A4 antibody [SM5058] (STJA0005058)

STJA0005058

GENERAL INFORMATION


| | |
|--------------------------|---|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-S100A4 is suitable for use in Immunohistochemistry and Immunofluorescence research applications. |
| Applications | IHC/IF |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | SM5058 |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen. |
| Dilution Range | IHC/IF (H) 1:200-1:500 |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Isotype | |
| Storage Instruction | Store at -20C for up to one year, and avoid repeated freeze-thaw cycles. |

TARGET INFORMATION

| | |
|-------------------------|---|
| Gene ID | 6275 |
| Gene Symbol | S100A4 |
| Uniprot ID | S10A4_HUMAN |
| Immunogen | Synthesized peptide derived from human S100A4 |
| Immunogen Region | |
| Specificity | |
| Immunogen | |
| Sequence | |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081