

## Anti-Cytokeratin 17 antibody [SM5097] (STJA0005097)

STJA0005097

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Cytokeratin 17 is suitable for use in Western Blot, Immunohistochemistry and Immunofluorescence research applications.
<b>Applications</b>	WB/IHC/IF
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	SM5097
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB (H) 1:500-1:2000 IHC/IF (H) 1:500-1:1000
<b>Formulation</b>	PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol.
<b>Isotype</b>	IgG2ak
<b>Storage</b>	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.
<b>Instruction</b>	

### TARGET INFORMATION

<b>Gene ID</b>	3872
<b>Gene Symbol</b>	KRT17
<b>Uniprot ID</b>	K1C17_HUMAN
<b>Immunogen</b>	KLH conjugated Synthetic peptide corresponding to Human Cytokeratin 17
<b>Immunogen Region</b>	
<b>Specificity</b>	
<b>Immunogen Sequence</b>	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081