

Anti-Aquaporin 3 antibody [SM5102] (STJA0005102)

STJA0005102

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Aquaporin 3 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Mouse
Reactivity	Mouse/Rat

PRODUCT PROPERTIES


Clonality	Monoclonal
Clone ID	SM5102
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	IHC/IF (M, R) 1:1000-1:2000
Formulation	PBS with 0.02% sodium azide, 100 µg/ml BSA and 50% glycerol.
Isotype	IgG
Storage Instruction	Store at -20°C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11828
Gene Symbol	Aqp3
Uniprot ID	AQP3_MOUSE
Immunogen	Recombinant protein corresponding to Mouse Aquaporin 3
Immunogen Region	
Specificity	
Immunogen Sequence	



Immunohistochemistry of paraffin embedded mouse kidney using AQP3 at dilution of 1:2000



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081