

Anti-ACTC1 antibody [SM5130] (STJA0005130)

STJA0005130

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-ACTC1 is suitable for use in Immunohistochemistry and Immunofluorescence research applications.
Applications	IHC/IF
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM5130
Concentration	
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution Range	IHC/IF (H) 1:200-1:500
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Isotype	
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	70 72 58 ACTC1 ACTG2 ACTC_HUMAN ACTH_HUMAN ACTS_HUMAN
Immunogen	Synthesized peptide derived from human Actin, Muscle Specific
Immunogen Region	
Specificity	
Immunogen Sequence	



Immunohistochemistry analysis of paraffin embedded human cardiac muscle using ACTC1 at dilution of 1:200