

Anti-CD41/CD61 antibody [SM5150] (STJA0005150)

STJA0005150

GENERAL INFORMATION


| | |
|--------------------------|--|
| Product Type | Primary antibodies |
| Short Description | Mouse monoclonal antibody anti-CD41/CD61 is suitable for use in Immunohistochemistry and Immunofluorescence research applications. |
| Applications | IHC/IF |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|--|
| Clonality | Monoclonal |
| Clone ID | SM5150 |
| Concentration | |
| Conjugation | Unconjugated |
| Purification | The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen. |
| Dilution Range | IHC/IF (H) 1:200-1:500 |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Isotype | |
| Storage Instruction | Store at -20C for up to one year, and avoid repeated freeze-thaw cycles. |

TARGET INFORMATION

| | |
|---------------------------|---|
| Gene ID | 3690 |
| Gene Symbol | ITGB3 |
| Uniprot ID | ITB3_HUMAN |
| Immunogen | Synthesized peptide derived from human CD61 |
| Immunogen Region | |
| Specificity | |
| Immunogen Sequence | |



Immunohistochemistry analysis of paraffin embedded human spleen using CD41/CD61 at dilution of 1:200

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081