

## Anti-Interferon gamma antibody [SM5198] (STJA0005198)

STJA0005198

### GENERAL INFORMATION


|                          |                                                                                                            |
|--------------------------|------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                         |
| <b>Short Description</b> | Mouse monoclonal antibody anti-Interferon gamma is suitable for use in Western Blot research applications. |
| <b>Applications</b>      | WB                                                                                                         |
| <b>Host/Source</b>       | Mouse                                                                                                      |
| <b>Reactivity</b>        | Human                                                                                                      |

### PRODUCT PROPERTIES

|                            |                                                                          |
|----------------------------|--------------------------------------------------------------------------|
| <b>Clonality</b>           | Monoclonal                                                               |
| <b>Clone ID</b>            | SM5198                                                                   |
| <b>Concentration</b>       |                                                                          |
| <b>Conjugation</b>         | Unconjugated                                                             |
| <b>Purification</b>        | Affinity purification                                                    |
| <b>Dilution Range</b>      | WB (H) 1:500-1:1000                                                      |
| <b>Formulation</b>         | PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol.           |
| <b>Isotype</b>             | IgG                                                                      |
| <b>Storage Instruction</b> | Store at -20C for up to one year, and avoid repeated freeze-thaw cycles. |

### TARGET INFORMATION

|                           |                                                                   |
|---------------------------|-------------------------------------------------------------------|
| <b>Gene ID</b>            | <a href="#">15978</a>                                             |
| <b>Gene Symbol</b>        | <a href="#">Ifng</a>                                              |
| <b>Uniprot ID</b>         | <a href="#">IFNG_MOUSE</a>                                        |
| <b>Immunogen</b>          | KLH conjugated Synthetic peptide corresponding to Mouse IFN gamma |
| <b>Immunogen Region</b>   |                                                                   |
| <b>Specificity</b>        |                                                                   |
| <b>Immunogen Sequence</b> |                                                                   |



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081