

Anti-FITC antibody [SM5205] (STJA0005205)

STJA0005205

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-FITC is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	FITC

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	SM5205
Concentration	
Conjugation	Unconjugated
Purification	None
Dilution Range	Elisa (Non) 1:5000-1:10000
Formulation	PBS with 0.02% sodium azide, 100 Mu g/ml BSA and 50% glycerol.
Isotype	IgG1k
Storage Instruction	Store at -20C for up to one year, and avoid repeated freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthesized peptide derived from FITC
Immunogen Region	
Specificity	
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081