

Anti-ABL1 antibody (C-Term) [M209] (STJA0005238)

STJA0005238

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-c-Abl (C-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M209
Concentration	
Conjugation	Unconjugated
Purification	Protein G Purified
Dilution Range	WB 1:1000
Formulation	PBS + 0.2% gelatin and 0.09% NaN ₃ .
Isotype	IgG1
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	25
Gene Symbol	ABL1
Uniprot ID	ABL1_HUMAN
Immunogen	Clone (M209) was generated from a recombinant protein corresponding to the C-terminal region of human c-Abl.
Immunogen Region	C-Term
Specificity	This antibody detects a 145 kDa* protein corresponding to c-Abl on SDS-PAGE immunoblots of human K-562 and HL-60 cells, and mouse ANN-1 cells.
Immunogen Sequence	