

Anti-ADIPOQ antibody (C-Term) [M576] (STJA0005247)

STJA0005247

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Adiponectin (C-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M576
Concentration	
Conjugation	Unconjugated
Purification	Protein G Purified
Dilution Range	WB 1:500
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	11450
Gene Symbol	Adipoq
Uniprot ID	ADIPO_MOUSE
Immunogen	Clone (M576) was generated from a recombinant protein that included amino acid residues within the C-terminal region of mouse Adiponectin. This sequence has high homology with similar regions in rat and human adiponectin.
Immunogen Region	C-Term
Specificity	This antibody detects a 30 kDa* protein on SDS-PAGE immunoblots of human and mouse serum, as well as detects full length human and mouse recombinant proteins.
Immunogen Sequence	