

Anti-Phospho-Pycard-Tyr-144 antibody (STJA0005271)

STJA0005271

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Asc-Tyr-144 is suitable for use in Western Blot and Immunocytochemistry research applications.
Applications	WB/ICC
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution	WB 1:1000
Range	ICC 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	66824
Gene Symbol	Pycard
Uniprot ID	ASC_MOUSE
Immunogen	Asc (Tyr-144) phospho-peptide (coupled to KLH) corresponding to amino acid residues surrounding Tyr-144 in mouse Asc. This peptide sequence is highly conserved in human and rat Asc.
Immunogen Region	
Specificity	This antibody was cross-adsorbed to unphosphorylated Asc (Tyr-144) peptide before affinity purification using phospho-Asc (Tyr-144) peptide. The antibody detects phosphorylated Asc (Tyr-144) in the inflammasome. In J774 macrophage cells primed with L
Immunogen Sequence	