

Anti-CDH1 antibody (C-Term) [M533] (STJA0005283)

STJA0005283

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-E-Cadherin (C-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M533
Concentration	
Conjugation	Unconjugated
Purification	Protein G Purified
Dilution Range	WB 1:250
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	12550
Gene Symbol	Cdh1
Uniprot ID	CADH1_MOUSE
Immunogen	Clone (M533) was generated from a mouse recombinant E-Cadherin protein containing amino acids in the C-terminal region. This sequence is highly conserved in human and rat E-cadherin.
Immunogen Region	C-Term
Specificity	The E-cadherin antibody detects a 120 kDa* protein corresponding to endogenous E-Cadherin expressed in human A431 cells.
Immunogen Sequence	