

Anti-CDH3 antibody (N-Term) [M596] (STJA0005291)

STJA0005291

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-P-Cadherin (N-Term) is suitable for use in Western Blot and Immunocytochemistry research applications.
Applications	WB/ICC
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M596
Concentration	
Conjugation	Unconjugated
Purification	Protein A Purified
Dilution Range	WB 1:1000 ICC 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1001
Gene Symbol	CDH3
Uniprot ID	CDH3_HUMAN
Immunogen	Clone (M596) was generated from a recombinant protein containing amino acids in the N-terminal region of human P-Cadherin. This sequence is highly conserved in mouse and rat P-cadherin.
Immunogen Region	N-Term
Specificity	The P-cadherin antibody detects a 120 kDa* protein corresponding to endogenous P-Cadherin expressed in human A431 cells.
Immunogen Sequence	