

Anti-CD244 antibody (Extracellular region) [M085] (STJA0005335)

STJA0005335

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-CD244/2B4/SLAMF4 M085 (Extracellular region) is suitable for use in Western Blot-Nat, Immunocytochemistry and Immunoprecipitation research applications.
Applications	WB-Nat/ICC/IP
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M085
Concentration	
Conjugation	Unconjugated
Purification	Protein G Purified
Dilution	WB 1:1000
Range	ICC 1:200 IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	51744
Gene Symbol	CD244
Uniprot ID	CD244_HUMAN
Immunogen	Clone (M085) was generated from a recombinant protein that included the extracellular region of human CD244 protein.
Immunogen Region	Extracellular region
Specificity	Clone M085 mouse monoclonal antibody detects a 110 kDa* protein on SDS-PAGE Native* immunoblots of human THP1 monocytes. The antibody detects only the native CD244 protein, and the aldehyde-fixed CD244 protein. The antibody works for native western b
Immunogen Sequence	