

Anti-CFL1 antibody (N-Term) (STJA0005345)

STJA0005345

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cofilin 1 (N-Term) is suitable for use in Western Blot and Immunocytochemistry research applications.
Applications	WB/ICC
Host/Source	Rabbit
Reactivity	Chicken/Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution	WB 1:1000
Range	ICC 1:50
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1072
Gene Symbol	CFL1
Uniprot ID	COF1_HUMAN
Immunogen	A synthetic peptide (coupled to KLH) corresponding to amino acid residues at the N-terminus of human Cofilin 1. This sequence has 100% homology with similar regions of rat and mouse Cofilin 1, and has 3 amino acid differences from human Cofilin 2.
Immunogen Region	N-Term
Specificity	The antibody detects a 19 kDa* protein corresponding to the molecular mass of Cofilin 1 on SDS-PAGE immunoblots of Jurkat cells.
Immunogen Sequence	