

Anti-FHOD1 antibody (Full Length) [M352] (STJA0005421)

STJA0005421

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-FHOD1 (Full Length) is suitable for use in Western Blot and Immunohistochemistry research applications.
Applications	WB/IHC
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M352
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution	WB 1:250
Range	IHC 1:50
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	29109
Gene Symbol	FHOD1
Uniprot ID	FHOD1_HUMAN
Immunogen	FHOD1 recombinant protein that included the full length human FHOD1 sequence. This sequence has 50% homology to human FHOD3, and low homology to other formin family members.
Immunogen Region	Full Length
Specificity	The antibody detects a 140 kDa* protein corresponding to the apparent molecular mass of FHOD1 on SDS-PAGE immunoblots of human K562 and Jurkat, and has weak reactivity to FHOD1 in mouse C2C12.
Immunogen Sequence	