

Anti-N/A antibody (C-Term) [M050] (STJA0005435)

STJA0005435

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-His Tag (C-Term) is suitable for use in Western Blot, Immunocytochemistry, Immunohistochemistry, Immunoprecipitation and Flow Cytometry research applications.
Applications	WB/ICC/IHC/IP/FC
Host/Source	Mouse
Reactivity	

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M050
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution	WB 1:5000
Range	IHC 1:100 ICC 1:100 IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Clone (M050) was generated from a recombinant protein containing a C-terminal His-tag.
Immunogen Region	C-Term
Specificity	Clone M050 detects recombinant proteins with C-terminal His-tag epitope. The antibody binds 6 and 10 His-tag epitopes on the C-terminal side of a protein, but does not bind N-terminal His-tag epitopes. The antibody works for western blot, immunopreci
Immunogen Sequence	