

Anti-PXN antibody [M107] (STJA0005512)

STJA0005512

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Paxillin is suitable for use in Western Blot, Immunocytochemistry and Immunoprecipitation research applications.
Applications	WB/ICC/IP
Host/Source	Mouse
Reactivity	Chicken/Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M107
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution	WB 1:1000
Range	ICC 1:100 IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5829
Gene Symbol	PXN
Uniprot ID	PAXI_HUMAN
Immunogen	Clone (M107) was generated from full-length recombinant human paxillin.
Immunogen Region	
Specificity	This antibody detects a 68kDa* protein corresponding to the molecular mass of paxillin on SDS-PAGE immunoblots of A431 cells. Similar results were seen in various cells and tissues from rat, mouse, and chicken origin.
Immunogen Sequence	