

Anti-PXN antibody (C-Term) (STJA0005513)

STJA0005513

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Paxillin (C-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	5829
Gene Symbol	PXN
Uniprot ID	PAXI_HUMAN
Immunogen	A synthetic peptide (coupled to KLH) corresponding to amino acid residues in the C-terminal region of human Paxillin. This peptide sequence is highly conserved in rat and mouse paxillins, and is found in all three paxillin isoforms (Alpha, Beta, and C-Term
Immunogen Region	C-Term
Specificity	This antibody detects a 68kDa* protein corresponding to the molecular mass of paxillin on SDS-PAGE immunoblots of A431 cells.
Immunogen Sequence	