

Anti-Phospho-N/A-pThr antibody (STJA0005527)

STJA0005527

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Phospho-Phosphothreonine-pThr is suitable for use in Western Blot, Immunocytochemistry and Immunoprecipitation research applications.
Applications	WB/ICC/IP
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution	WB 1:1000
Range	ICC 1:250 IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Phosphothreonine synthetic peptides coupled to carrier protein.
Immunogen Region	
Specificity	This antibody was cross-adsorbed to unphosphorylated peptide then affinity purified using phospho-threonine peptide (without carrier). The antibody detects many serine or threonine phosphorylated proteins by western blot and ELISA.
Immunogen Sequence	