

Anti-PLXND1 antibody (C-Term) (STJA0005544)

STJA0005544

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Plexin D1 (C-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	23129
Gene Symbol	PLXND1
Uniprot ID	PLXD1_HUMAN
Immunogen	Plexin D1 synthetic peptide (coupled to carrier protein) corresponding to amino acids from the C-terminal region of human Plexin D1. This sequence is well conserved in rat and mouse Plexin D1, and has some homology to other Plexin family members.
Immunogen Region	C-Term
Specificity	The antibody detects 215 and 195 kDa bands corresponding to the apparent molecular mass of Plexin D1 variants in human umbilical vein endothelial cells.
Immunogen Sequence	