

Anti-SRC antibody (N-Term) [M259] (STJA0005601)

STJA0005601

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-c-Src (N-Term) is suitable for use in Western Blot and Immunoprecipitation research applications.
Applications	WB/IP
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	M259
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution	WB 1:1000
Range	IP 1:100
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG1
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	6714
Gene Symbol	SRC
Uniprot ID	SRC_HUMAN
Immunogen	Clone (M259) was generated from a recombinant protein containing amino acids in the N-terminal region of human c-Src. This sequence has high homology to similar regions in rat and mouse c-Src.
Immunogen Region	N-Term
Specificity	This antibody detects a 60 kDa* protein corresponding to c-Src on SDS-PAGE immunoblots of mouse SYF cells transformed with c-Src and human A431 cells.
Immunogen Sequence	