

Anti-FLT4 antibody (N-Term) (STJA0005642)

STJA0005642

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-VEGFR-3 (N-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:1000
Formulation	PBS + 1 mg/ml BSA, 0.05% NaN ₃ and 50% glycerol
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2324
Gene Symbol	FLT4
Uniprot ID	VEGFR3_HUMAN
Immunogen	A synthetic peptide (coupled to carrier protein) corresponding to amino acids from the N-terminal region of human VEGFR-3. This sequence has two amino acid differences from rat VEGFR-3, and five amino acid differences from mouse VEGFR-3.
Immunogen Region	N-Term
Specificity	The antibody detects a 170 kDa* band corresponding to VEGFR-3 in western blots of human endothelial and K-562 cells, and shows strong reactivity to recombinant human VEGFR-3, but not VEGFR-1 or VEGFR-2.
Immunogen Sequence	