

Anti-Adenovirus 3 antibody (STJA0005664)

STJA0005664

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Adenovirus 3 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human Adenovirus 3

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN3
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Synthetic peptide corresponding to amino acids from the fiber knob protein of human Ad 3.
Immunogen Region	
Specificity	Reacts with an antigenic epitope of the fiber knob of Ad Serotype 3. The fiber knob is involved in the initial interaction of adenovirus with cellular membrane receptors.
Immunogen Sequence	