

Anti-PTCH1 antibody (STJA0005679)

STJA0005679

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-PTCH1 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Mouse/Chicken/Human/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	19206
Gene Symbol	Ptch1
Uniprot ID	PTC1_MOUSE
Immunogen	Synthetic peptide corresponding to amino acids from mouse PTCH1 protein
Immunogen Region	
Specificity	Reacts with mouse PTCH protein. Expected to cross react with rat, chicken and human due to 100% sequence identity.
Immunogen Sequence	