

Anti-CD28 antibody (Ig-V-like Region) (STJA0005685)

STJA0005685

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-CD28 (Ig-V-like Region) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	940
Gene Symbol	CD28
Uniprot ID	CD28_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids from the Ig-V-like region human CD28 protein
Immunogen Region	Ig-V-like Region
Specificity	Reacts with the Ig-V-like region of human CD28 protein.
Immunogen Sequence	