

Anti-IKBKG antibody (N-Term) (STJA0005690)

STJA0005690

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-IKK gamma (N-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	8517
Gene Symbol	IKBKG
Uniprot ID	NEMO_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids from the N-terminal region of human IKK gamma protein
Immunogen Region	N-Term
Specificity	Reacts with the N-terminal region of human IKK gamma protein. Recognizes both human isoforms.
Immunogen Sequence	