

**Anti-ELP1 antibody (N-Term) (STJA0005692)**

STJA0005692

**GENERAL INFORMATION**

<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Rabbit polyclonal antibody anti-IKAP (N-Term) is suitable for use in Western Blot research applications.
<b>Applications</b>	WB
<b>Host/Source</b>	Rabbit
<b>Reactivity</b>	Human

**PRODUCT PROPERTIES**

<b>Clonality</b>	Polyclonal
<b>Clone ID</b>	
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Antigen Affinity Purified
<b>Dilution Range</b>	WB 1:500-1:1000
<b>Formulation</b>	PBS + 0.02% NaN <sub>3</sub>
<b>Isotype</b>	IgG
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

**TARGET INFORMATION**

<b>Gene ID</b>	<a href="#">8518</a>
<b>Gene Symbol</b>	<a href="#">ELP1</a>
<b>Uniprot ID</b>	<a href="#">ELP1_HUMAN</a>
<b>Immunogen</b>	Synthetic peptide corresponding to amino acids from the N-terminal region of human IKAP protein
<b>Immunogen Region</b>	N-Term
<b>Specificity</b>	Reacts with the N-terminus of human IKAP protein. Sequence is 92% identical to rat and 85% identical to mouse.
<b>Immunogen Sequence</b>	