

Anti-SHANK1 antibody (N-Term) (STJA0005704)

STJA0005704

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-SHANK1a (N-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	
Conjugation	Unconjugated
Purification	Antigen Affinity Purified
Dilution Range	WB 1:500-1:1000
Formulation	PBS + 0.02% NaN ₃
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	50944
Gene Symbol	SHANK1
Uniprot ID	SHAN1_HUMAN
Immunogen	Synthetic peptide corresponding to amino acids from the N-terminal region of human Shank1a protein
Immunogen Region	N-Term
Specificity	Reacts with the N-terminal region of human, mouse and rat Shank1a protein.
Immunogen Sequence	